

SCA 88

strong alkaline laboratory detergent

SCA LABORATORY DETERGENT SERIES

a product of



SCA

SCA PACIFIC Pte Ltd

SCA 88 STRONG ALKALINE LABORATORY DETERGENTS



SCA 88 Strong Alkaline Laboratory Detergent is specially formulated for use in the different type of laboratory glassware washing machine and different type of ultrasonic bath. The Strong Alkaline Laboratory Detergent contains non-ionic, an-ionic, and phosphate **OR WITHOUT** phosphate surfactants to give good emulsifying and dispersing action as well as with active chlorine as disinfectant component provide bactericide, fungicide, turboculocide and virus inactivating action.

Direction of used: 1:10 dilution with deionizer water

Packing Size: 1 liter and 5 liters

SCA Pacific Pte Ltd

Address: Blk 108 Ang Mo Kio Ave 4, #01-88
Singapore 560108

Tel: +65 6298 9989

Fax: +65 6291 9989

E-mail: scapacific@singnet.com.sg

Website: www.scapacific.com.sg

Contact Our Authorized Distributor

www.scapacific.com.sg

In line with our policy of continuous improvement, SCA Pacific Pte Ltd reserves the right to revise published specifications, technical information and details of the products without prior notice.